

# Comparative Study of DNA Isolated from Whole Blood Using Norgen's Blood Genomic DNA Isolation Mini Kit and Competitors' Kits

M.EL-MOGY, PH.D<sup>1</sup>, M. SIMKIN<sup>1</sup>, Y. HAJ-AHMAD, PH.D<sup>1,2</sup>

<sup>1</sup>NORGEN BIOTEK CORPORATION, THOROLD, ONTARIO, CANADA

<sup>2</sup>CENTRE FOR BIOTECHNOLOGY. BROCK UNIVERSITY, ST. CATHARINES, ONTARIO, CANADA

## INTRODUCTION

Blood is an excellent sample for diagnostic purposes. It is rich in protein and nucleic acid-based biomarkers that can be used to detect a panel of diseases long before any physical symptoms become apparent. Studies have shown that cancer biomarkers, such as differentially methylated target genes, are detectable in blood samples and have proven to be extremely sensitive and specific for given cancer types<sup>1</sup>. These biomarkers can be cell associated, such as those found in leukocytes<sup>2</sup>, or cell-free, such as those found in plasma or serum samples. For this reason, blood is often the sample of choice for biomarker or diagnostic research.

Investigators utilizing blood in their research have unique needs, based on their downstream applications and sample volume. The method used for blood DNA isolation can heavily influence the results of research-based or diagnostic tests associated with blood. Investigators must ensure that their blood DNA isolation method is efficient in isolating high amounts of DNA with high purity and compatibility with sensitive downstream applications. A robust blood DNA isolation kit eliminates sample processing biases, and increases data reproducibility.

The purpose of this study is to compare four commercially available blood DNA isolation methods: Norgen's Blood Genomic DNA Isolation Mini Kit (Cat# 46300), Competitor A's DNA Blood Mini Kit), Competitor B's Blood kit and Competitor B's Dx Blood, for their ability to isolate high quality and high quantities of genomic DNA from increasing volumes of blood.



## METHODS AND MATERIALS

### Sample collection

Ten milliliters of blood was collected in Ethylenediaminetetraacetic acid (EDTA) tubes, from one healthy individual, by a trained professional. The sample was frozen at -70°C until processed.

### Blood DNA extraction

DNA was extracted from 200 µL of the thawed blood sample using Norgen's Blood Genomic DNA Isolation Mini Kit (Cat# 46300), Competitor A's DNA Blood Mini Kit, Competitor B's Blood kit and Competitor B's Dx Blood, as per the manufacturers' instructions.

### Spectrophotometry

The purified RNA was then used as the template in both real-time PCR and real-time RT-PCR reactions. In both sets of reactions, primers that recognize both 16S rRNA and 16S rDNA from prokaryotes were used for amplification at a final concentration of 0.25 µM, and 2 µL of RNA elution was added as template. SYBR Green was used for detection in both sets of reactions. Recombinant Taq Polymerase was incorporated into both the PCR and the RT-PCR reactions and SuperScript III reverse transcriptase was incorporated into only the one step real-time RT-PCR reactions. One step RT-qPCR; 50°C for 30 minutes for reverse transcription, 95°C for 3 minutes for denaturation, and 45 cycles of 95°C for 15 seconds for denaturation, 60°C for 30 seconds for annealing and 72°C for 45 seconds for extension. Meltingcurve analysis was also performed. The reaction was run on an iCycler iQ real-time system (Bio-Rad).


### Real-Time PCR

The purified DNA was then used as the template in a realtime TaqMan® PCR reaction. Briefly, 1, 3, 6 & 9 µL of isolated DNA was added to 20 µL of real-time PCR reaction mixture containing 10 µL of Norgen's 2X PCR Mastermix (Cat# 28007), 0.4 µL from a 25mM GAPDH primer pair mix, 0.2 µL from a 25mM TaqMan® probe, and nuclease-free water. The PCR samples were amplified under the real-time program; 95°C for 3 minutes for an initial denaturation, 40 cycles of 95°C for 15 seconds for denaturation, 60°C for annealing and extension. The reaction was run on an iCycler iQ Realtime System (Bio-Rad).

## RESULTS AND DISCUSSION

Blood is an excellent resource for research and diagnostic purposes. Systemic diseases can be detected through a simple blood test, and changes in DNA (both cellular and cell-free) can be detected from blood, giving rise to its high potential for screening for a panel of diseases.


The key to the success of any study utilizing blood is a reliable blood DNA isolation method that provides high quality DNA. In this study, DNA was isolated from 200  $\mu$ L of blood using Norgen's kit, Competitor A and two Competitor B kits.



**Figure 1. Resolution of DNA isolated using Norgen's Blood Genomic DNA Isolation Kit (Norgen) and Competitor A's DNA Blood Mini Kit (A), Competitor B's Blood kit (B) and Competitor B's Dx Blood (B Dx). Fifteen microliters of 200  $\mu$ L elutions were run on 1X TAE 1.0% agarose gel. Marker= Norgen's UltraRanger DNA Ladder.**

Fifteen microliters of each 200  $\mu$ L elution was then run on a 1X TAE 1.0% agarose gel to visually inspect the isolated gDNA (**Figure 1**). All kits were able to isolate intact genomic DNA of > 24kbp in size, with higher yields from Norgen's Blood Genomic DNA Isolation Mini kit and Competitor B's Dx Blood kit.

In order to determine the yield and purity of the isolated DNA, samples were measured using a cuvette-based spectrophotometry method (UltraSpec 2100 Pro; Fisher Scientific) (**Figure 2 and Figure 3**). Norgen samples were found to have higher DNA yields and OD260/280 ratios than the other competitors with an average recovery of 4  $\mu$ g of DNA and OD260/280 ratio of 2.



**Figure 2.** The difference in yield of blood DNA isolated using Norgen's Blood Genomic DNA Isolation Mini Kit (Norgen) and Competitor A's DNA Blood Mini Kit (A), Competitor B's Blood kit (B) and Competitor B's Dx Blood (B Dx). Fifty microliters of each sample was diluted in 450  $\mu$ L of nuclease-free water, and DNA concentrations were measured using the UltraSpec 2100 Pro (Fisher Scientific).

DNA quality was determined through the use of a TaqMan® Real-Time PCR method. In order to assess sample inhibition, 1, 3, 6 & 9 µL of each sample was used in the reaction. The Ct values were then graphed (Figure 4).

Successful amplification was obtained from all samples for all the input volumes tested. Norgen's kit showed the lowest average Ct value at 9 µL template input compared to the other kits, indicating the higher purity of the isolated DNA and absence of PCR inhibitors.




Figure 3. The difference in OD260/280 of blood DNA isolated using Norgen's Blood Genomic DNA Isolation Mini Kit (Norgen) and Competitor A's DNA Blood Mini Kit (A), Macherey-Nagel's Competitor B's Blood kit (B) and Competitor B's Dx Blood (B Dx). Fifty microliters of each sample was diluted in 450 µL of nuclease-free water, and DNA concentrations were measured using the UltraSpec 2100 Pro (Fisher Scientific).




Figure 4. The difference in Ct values between DNA isolated using Norgen’s Blood Genomic DNA Isolation Mini Kit (Norgen), Competitor A’s DNA Blood Mini Kit (A), Competitor B’s Blood kit (B) and Competitor B’s Dx Blood (B Dx), in a Taqman® qPCR reaction. Different template volumes (1, 3, 6 & 9  $\mu$ L) microliters of each elution were used in a 20 $\mu$ L qPCR reaction involving GAPDH primers.

## CONCLUSION

From the data presented in this report, the following can be concluded:

1. Norgen’s Blood Genomic DNA Isolation Mini Kit isolates the highest yields of DNA, when compared to the leading competitors’ kits.
2. Norgen’s Blood Genomic DNA Isolation Mini Kit outperforms competitors’ kits in sensitive downstream applications, such as qPCR.
3. Norgen samples are free from inhibitors, as 9  $\mu$ L of elution was used in each PCR reaction, with no delay in Ct values when compared to competitors’ kits.

## REFERENCES

1. Warren JD, Xiong W, Bunker AM, Vaughn CP, Furtado LV, Roberts WL, et al. 2011. Septin 9 methylated DNA is a sensitive and specific blood test for colorectal cancer. *BMC Med*; 9 (133).
2. Koestler DC, Marsit CJ, Christensen BC, Accomando W, Langevin SM, Houseman EA, et al. 2012. Peripheral blood immune cell methylation profiles are associated with nonhematopoietic cancers. *Cancer Epidemiol Biomarkers Prev*; 21(8):1293-302.



*Norgen is an ISO certified company,  
indicating our commitment to quality*



## Order Today!

### Three ways to Order

Visit [norgenbiotech.com](http://norgenbiotech.com)

Call at **1-866-NORGENB** or **905-227-8848**

Email at [orders@norgenbiotech.com](mailto:orders@norgenbiotech.com)

## RELATED PRODUCTS

Description	Size	Cat. #
Blood DNA Isolation Kits	50 Preps	46300



For more information  
**SCAN HERE**