

FBS Exosome Depletion Kits

CAT. 61100, 61400, 61200, 61300

QUICK AND EASY PROTOCOL FOR THE DEPLETION OF BOVINE EXOSOMES FROM FBS PRIOR TO USING IT AS A GROWTH SUPPLEMENT IN YOUR CULTURE MEDIUM




- ✓ Efficient depletion of bovine exosomes from Fetal Bovine Serum
- ✓ Deplete exosome-sized vesicles from versatile FBS volumes
- ✓ No protease treatment required, no time-consuming ultracentrifugation, no precipitation reagents required, no overnight incubation required
- ✓ Exosome depletion confirmed by reduction of bovine miRNAs below detectable levels
- ✓ The depleted FBS provides the same cellular growth rates as the standard FBS
- ✓ Purification is based on Norgen's proprietary Silicon Carbide resin




Advantages You Will Bring to Your Lab

A



B



RT-PCR ANALYSIS OF BOVINE MIRNA LEVELS SHOWS DEPLETION OF EXOSOMES USING NORGEN'S SLURRY BASED FBS EXOSOME DEPLETION KIT


Figure 1. Exosome-depleted FBS with Norgen's FBS Exosome Depletion Kits (Slurry Format) has undetectable Bovine miRNA levels. Norgen's FBS Exosome Depletion Kit I (Slurry Format) (Cat# 61100) was used to deplete bovine exosomes from 5mL FBS. Total RNA/miRNA including exosomal RNA was purified from the depleted FBS, non-depleted FBS and a commercially available ready – to – go depleted FBS. Five different bovine microRNAs were assessed by RT-qPCR (miR-26a, miR-30a, miR-92a, miR-23a and miR-122). Three out of the five tested miRNA (miR-26a, miR-92a and miR-23a) didn't show any amplification in the FBS depleted using Norgen's FBS Exosome Depletion Kit I (Slurry Format) whereas the other two miRNAs (miR-30a and miR-122) showed very late Ct. values which appeared to be a primer dimer according to the melt curve.

FBS EXOSOME DEPLETION DOESN'T AFFECT CELLULAR MORPHOLOGY AND GROWTH


Figure 2. Growth rates of HeLa cells in media containing Exosome-depleted FBS.

Growth rates of HeLa cells in media containing Exosome-depleted FBS using Norgen's FBS Exosome Depletion Kits (Slurry Format) was compared to that in media containing standard FBS. Simply, HeLa cells were seeded in DMEM with either 10% Exosome-depleted FBS using Norgen's Kits or 10% standard FBS and then cultured under standard conditions at 37°C with 5% CO₂ for 3 days. The cells were imaged using Moticam 480 to observe cellular morphology and growth rate. Similar growth and identical cellular morphology were detected for both the Exosome-depleted FBS using Norgen's FBS Exosome Depletion Kits and the standard FBS.


Standard FBS




Depleted FBS Norgen



A



B



RT-PCR ANALYSIS OF BOVINE MIRNA LEVELS SHOWS DEPLETION OF EXOSOMES USING NORGEN'S COLUMN BASED FBS EXOSOME DEPLETION KIT

Figure 3. Exosome-depleted FBS with Norgen's FBS Exosome Depletion Kits (Column Format) has undetectable Bovine miRNA levels. Norgen's FBS Exosome Depletion Kit I (Column Format) (Cat# 61200) was used to deplete bovine exosomes from 5mL FBS. Total RNA/miRNA including exosomal RNA was purified from the depleted FBS, non-depleted FBS and a commercially available ready – to – go depleted FBS. Five different bovine microRNAs were assessed by RT-qPCR (miR-26a, miR-30a, miR-92a, miR-23a and miR-122). Three out of the five tested miRNA (miR-26a, miR-92a and miR-23a) didn't show any amplification in the FBS depleted using Norgen's FBS Exosome Depletion Kit I (Column Format) whereas the other two miRNAs (miR-30a and miR-122) showed very late Ct. values which appeared to be a primer dimer according to the melt curve.

TECHNICAL SPECIFICATIONS

Description	Specifications
Sample Type	Fetal Bovine Serum
Sample Volume Range	Up to 140 mL (FBS Exosome Depletion Kit I (Slurry Format)) Up to 280 mL (FBS Exosome Depletion Kit II (Slurry Format)) Up to 120 mL (FBS Exosome Depletion Kit I (Column Format)) Up to 240 mL (FBS Exosome Depletion Kit II (Column Format))
Depletion	Deplete exosome-sized vesicle
Bovine miRNA	No detectable bovine miRNA
Time to Complete 6 Purifications	40 minutes

SELECT PUBLICATIONS

Publication Title	Authors	Journal	Year
Excessive activation of IL-33/ST2 in cancer-associated fibroblasts promotes invasion and metastasis in ovarian cancer	Caixia Feng, Li Kou, Panyue Yin, Yuan Jing	Oncology Letters	2022
Tridimensional cell culture of dermal fibroblasts promotes exosome-mediated secretion of extracellular matrix proteins	Vincent Clément, Vincent Roy, Bastien Paré, Cassandra R. Goulet, Lydia Touzel Deschênes, François Berthod, Stéphane Bolduc & François Gros-Louis	Scientific Reports	2022

Ordering Information

Description	Preps	Cat. #
FBS Exosome Depletion Kit I	6 Preps	61100
FBS Exosome Depletion Kit II	12 Preps	61400
FBS Exosome Depletion Kit Column Format	6 Preps	61200
FBS Exosome Depletion Kit Column Format	12 Preps	61300

Order Today!

Order Norgen's FBS Exosome Depletion Kit I today to help improve your research performance and workflow.

Three ways to order

Visit norgenbiotek.com

Call at **1-866-NORGENB** or **905-227-8848**

Email at orders@norgenbiotek.com

Enhance your lab with the complete workflow solutions from Norgen.



COLLECTION & PRESERVATION



EXTRACTION



APPLICATION

Related Products	Prep Size	Cat. #
Exosomal RNA Isolation Kit	Variable	58000
Norgen Next Generation Sequencing Services	Call 1-866-NORGENB or Visit norgenbiotek.com	



For more information
SCAN HERE



norgenbiotek.com