

RNA Clean-Up and Concentration Micro-Elute Kit - Supplementary Protocol

Product # 61000

NGS Library Clean-Up and Concentration (removal of adapters, salts and enzymes)

Notes Prior to Use

- All centrifugation steps are carried out in a benchtop microcentrifuge at 14,000 x g (~14,000 RPM) except where noted. All centrifugation steps are performed at room temperature.
- A variable speed centrifuge should be used for maximum kit performance. If a variable speed centrifuge is not available a fixed speed centrifuge can be used, however reduced yields may be observed.
- Ensure that all solutions are at room temperature prior to use.
- Prepare a working concentration of the Wash Solution A by adding 90 mL of 96 - 100% ethanol (provided by the user) to the supplied bottle(s) containing the concentrated Wash Solution A. This will give a final volume of 128 mL. The label on the bottle has a box that may be checked to indicate that the ethanol has been added.
- It is important to work quickly during this procedure.

1. Sample Preparation

- a. Adjust the volume of the sample to 100 μ L by adding RNase-free water.
- b. Add 300 μ L of **Buffer RL** to the sample. Mix by vortexing
- c. Add 50 μ L of 96 – 100% ethanol (provided by the user) to the mixture from **Step 1b**. Mix by vortexing for 10 seconds.

2. Column Activation and Sample Binding to Column

- a. Assemble a column with one of the provided collection tubes.
- b. Apply 500 μ L of **Column Activation Solution** onto the column and centrifuge for 1 minute at 14,000 x g (~14,000 RPM).
- c. Discard the flowthrough. Reassemble the spin column with its collection tube.
- d. Apply the sample with the ethanol (from **Step 1c**) onto the **activated column** and centrifuge for 1 minute at $\geq 3,500$ x g (6,000 RPM).
- e. Discard the flowthrough. Reassemble the spin column with its collection tube.

3. Column Wash

- a. Apply 600 μ L of **Wash Solution A** to the column and centrifuge for 1 minute at 14,000 x g (~14,000 RPM).

Note: Ensure the entire Wash Solution A has passed through into the collection tube by inspecting the column. If the entire wash volume has not passed, spin for an additional minute.

- b. Discard the flowthrough and reassemble the spin column with its collection tube.
- c. Repeat steps **3a** and **3b** to wash the column a second time.
- d. Discard the flowthrough and the collection tube and reassemble the spin column with a new collection tube.
- e. Spin the column for 2 minutes at 14,000 x g (~14,000 RPM) in order to thoroughly dry the resin. Discard the collection tube.

4. DNA Elution

- a. Place the column into a fresh 1.7 mL Elution tube provided with the kit.
- b. Add 8 – 15 μ L of **Elution Solution A** to the column.

Note: For maximum concentrations, the 8 μ L elution volume may be used. For maximum recovery, the 15 μ L elution volume is recommended

- c. Centrifuge for 1 minute at **200 x g (~2,000 RPM)**, followed by 2 minutes at **14,000 x g (~14,000 RPM)** Note the volume eluted from the column. If the entire volume has not been eluted, spin the column at 14,000 x g (~14,000 RPM) for 1 additional minute.

Note: For maximum recovery, it is recommended that a second elution be performed into a separate microcentrifuge tube (Repeat **Steps 4b and 4c**).

5. Storage of DNA

The purified sample may be stored at -20°C .

Technical Support

Contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362.

Technical support can also be obtained from our website (www.norgenbiotek.com) or through email at techsupport@norgenbiotek.com.

3430 Schmon Parkway, Thorold, ON Canada L2V 4Y6
Phone: (905) 227-8848
Fax: (905) 227-1061
Toll Free in North America: 1-866-667-4362