

## Genomic DNA Isolation from 1 $\mu$ L – 100 $\mu$ L of Whole Blood using Norgen's Blood Genomic DNA Isolation Micro Kit

A. Tandon, PhD<sup>1</sup>, L. Graziano<sup>2</sup>, M. El-Mogy, PhD<sup>1</sup>, Y. Haj-Ahmad, Ph.D<sup>1,3</sup>

<sup>1</sup>Norgen Biotek Corporation, Thorold, Ontario, Canada

<sup>2</sup>Department of Biomedical Sciences, University of Guelph, Guelph, Ontario, Canada.

<sup>3</sup>Centre for Biotechnology, Brock University, St. Catharines, Ontario, Canada

### INTRODUCTION

Blood is an excellent sample for diagnostic purposes. It is rich in protein and nucleic acid-based biomarkers that can be used to detect a panel of diseases long before any physical symptoms become apparent. Studies have shown that cancer biomarkers, such as differentially methylated target genes, are detectable in blood samples and have proven to be extremely sensitive and specific for given cancer types<sup>1</sup>. These biomarkers can be cell associated, such as those found in leukocytes<sup>2</sup>, or cell-free, such as those found in plasma or serum samples. For this reason, blood is often the sample of choice for biomarker or diagnostic research.

Investigators utilizing blood in their research have unique needs, based on their downstream applications and sample volume. For some studies, less blood may be processed in order to preserve precious samples. On the other hand, some studies require a high yield of DNA from their samples, and thus require higher volumes of sample to be processed. The method used for blood DNA isolation can heavily influence the results of research-based or diagnostic tests associated with blood. Investigators must ensure that their blood DNA isolation method is flexible, i.e. a linear increase in blood volume being processed leads to a linear increase in DNA concentration and yield. A robust blood DNA isolation kit eliminates sample processing biases, and increases data reproducibility.

The purpose of this study is to test Norgen's Blood Genomic DNA Isolation Micro Kit (Cat# 52100) over a range of different blood input volumes (1  $\mu$ L - 100  $\mu$ L) for its ability to isolate high quality and high quantities of genomic DNA.

### MATERIALS AND METHODS

#### Sample collection

Blood was collected in Ethylenediaminetetraacetic acid (EDTA) tubes, from one healthy individual, by a trained professional. The sample was frozen at -70°C until processed.

#### Blood DNA extraction

DNA was extracted from the thawed blood samples using Norgen's Blood Genomic DNA Isolation Micro Kit (Cat# 52100) as per the manufacturer's instruction. Briefly, Proteinase K was added to a microcentrifuge tube, followed by 1, 10, 25, 50 or 100  $\mu$ L of blood. The volume of each was raised to 100  $\mu$ L using nuclease-free water.

Lysis Solution was then added, and samples were vortexed and incubated at 55°C for 10 minutes. Next, ethanol was added to each sample, and samples were bound, washed and eluted as per the manufacturer's protocol.

#### Spectrophotometry

Blood DNA quantity was measured using the UltraSpec 2100 Pro (Fisher Scientific). Fifty microliters of each DNA elution was diluted with 450  $\mu$ L of nuclease-free water, and OD measurements were taken using the cuvette-based spectrophotometry method.


#### Real-Time PCR

The purified DNA was then used as the template in a real-time TaqMan® PCR reaction. Briefly, 5  $\mu$ L of isolated DNA was added to 20  $\mu$ L of real-time PCR reaction mixture containing 10  $\mu$ L of Norgen's 2X PCR Mastermix (Cat# 28007), 0.4  $\mu$ L of 25  $\mu$ M GAPDH primer pair mix, 0.2  $\mu$ L of 25  $\mu$ M TaqMan® probe, and nuclease-free water. The PCR samples were amplified under the real-time program; 95°C for 3 minutes for an initial denaturation, 40 cycles of 95°C for 15 seconds for denaturation, 60°C for annealing and extension. The reaction was run on an iCycler iQ Realtime System (Bio-Rad).

### RESULTS AND DISCUSSION


Blood is an excellent resource for research and diagnostic purposes. Systemic diseases can be detected through a simple blood test, and changes in DNA (both cellular and cell-free) can be detected from blood, giving rise to its high potential for screening for a panel of diseases.

The key to the success of any study utilizing blood is a reliable blood DNA isolation method that performs optimally from a range of sample input volumes. One can determine the flexibility of a kit by increasing the sample input volume in order to see the linearity of the increase in DNA yield. In this study, DNA was isolated from 1, 10, 25, 50 and 100  $\mu$ L of blood using Norgen's kit. Fifteen microliters of each 100  $\mu$ L elution was then run on a 1X TAE 1.0% agarose gel to visually inspect the isolated genomic DNA (**Figure 1**). It was found that the yield was proportional to the blood input volume.




**Figure 1. Resolution of genomic DNA from different input volumes of blood that were collected in EDTA tubes and isolated using Norgen's Blood Genomic DNA Isolation Micro Kit.** Fifteen microliters of 100  $\mu$ L elutions were run on 1X TAE 1.0% agarose gel. Marker = Norgen's UltraRanger DNA Ladder.

In order to determine the yield and purity of DNA isolated from each volume of blood using Norgen's kit, samples were measured using a cuvette-based spectrophotometry method (Figure 2). The quantification did not agree with the gel, particularly at low input volumes, since the yield was lower than the detection limit of the device, however DNA was still visible on the gel. On the other hand, OD260/280 ratio was found to be > 1.7 from all samples (Figure 3).




**Figure 2. Genomic DNA yield isolated from different input volumes of blood using Norgen's Blood Genomic DNA Isolation Micro Kit.** Fifty microliters of each sample was diluted in 450 µL of nuclease-free water, and DNA concentrations were measured using the UltraSpec 2100 Pro (Fisher Scientific).



**Figure 3. OD260/280 ratios of DNA isolated from different input volumes of blood processed with Norgen's Blood Genomic DNA Isolation Micro Kit.** Fifty microliters of each sample was diluted in 450 µL of nuclease-free water, and DNA concentrations were measured using the UltraSpec 2100 Pro (Fisher Scientific).

DNA quality was determined through the use of a TaqMan® Real-Time PCR method. Five microliters of each sample was used in the reaction, and the Ct values were then graphed (Figure 4). Successful PCR amplification was observed for all samples with an almost linear decrease in Ct with increasing blood input volume, which indicates that the kit has a linear dynamic range for blood input ranges of 1 – 100 µL.



**Figure 4. The Ct values obtained from a Taqman® qPCR reaction performed on DNA isolated from different input volumes of blood using Norgen's Blood Genomic DNA Isolation Micro Kit.** Five microliters of each elution was used in a 20 µL qPCR reaction involving GAPDH primers.

## CONCLUSIONS

From the data presented in this report, the following can be concluded:

- 1- Isolated genomic DNA has good purity and integrity with positive PCR amplification from all samples over the blood input range tested (1 µL - 100 µL).
- 2- Norgen's Blood Genomic DNA Isolation Micro Kit displayed a good performance with a linear dynamic blood input range of 1 – 100 µL.

## REFERENCES

1. Warren JD, Xiong W, Bunker AM, Vaughn CP, Furtado LV, Roberts WL, et al. 2011. Septin 9 methylated DNA is a sensitive and specific blood test for colorectal cancer. *BMC Med*; 9 (133).
2. Koestler DC, Marsit CJ, Christensen BC, Accomando W, Langevin SM, Houseman EA, et al. 2012. Peripheral blood immune cell methylation profiles are associated with nonhematopoietic cancers. *Cancer Epidemiol Biomarkers Prev*; 21(8):1293-302.