

Total RNA Purification Plus Micro Kit

Norgen's Total RNA Purification Plus Micro Kit provides a rapid and sensitive method for the isolation and purification of total RNA from small input amounts of cultured animal cells, tissue samples, and microdissected samples including laser-capture microdissection (LCM). The kit purifies all sizes of RNA, from large mRNA and ribosomal RNA down to microRNA (miRNA) and small interfering RNA (siRNA). Genomic DNA is removed from the sample using a Genomic DNA Removal column, and the RNA is preferentially purified from other cellular components such as proteins, without the use of phenol or chloroform. The purified RNA is of the highest integrity, and can be used in a number of downstream applications including real time PCR, reverse transcription PCR,



Northern blotting, RNase protection and primer extension, and expression array assays.

Kit Specifications		
Column Binding Capacity	35 µg	Maximum Amount of Starting Material:
Maximum Column Loading Volume	650 µL	
Minimum Elution Volume	20 µL	Animal Tissues* 3 mg
Size of RNA Purified	All sizes	Laser-Captured Microdissection (LCM) Up to 5 x 10 ⁵ cells
Average RNA Yield:		Time to Complete 10 Purifications 22 minutes
HeLa Cells (1 x 10 ⁵ cells)	1.5 µg	

* for fibrous tissue, an additional Proteinase K treatment is required

Total RNA Purification Plus Micro Kit Benefits

Small elution volumes	The special design of the micro spin-column allows a small elution volume of as little as 20 µL.
Isolate a diversity of RNA species	All RNA species can be isolated, from large mRNA and ribosomal RNA down to microRNA (miRNA) and small interfering RNA (siRNA), without compromising total yield.
Fast and easy processing	Rapid spin-column format allows for the processing of 10 samples in 20 minutes.
No phenol:chloroform extractions	Total RNA is isolated without the use of harmful chemicals such as phenol or chloroform.
Ideal for small samples	Isolate total RNA from small input amounts of cultured animal cells, tissue samples, and LCM samples. Total RNA has been isolated and detected from as little as a single animal cell.
Recovered RNA is suitable for downstream applications	Purified RNA can be used in a number of downstream applications including real-time PCR, reverse transcription PCR, Northern blotting, RNase protection and primer extension, and expression array analysis requiring the use of intact RNA.

Total RNA Purification Plus Micro Kit




Figure 1. Great Isolation Sensitivity. Norgen's Total RNA Purification Plus Micro Kit allows sensitive RNA extraction from as little as a single cell. Total RNA was extracted from a decreasing number of HeLa cells followed by RT-qPCR to detect the human S14 transcript in the isolated RNA. PCR product of S14 was detected from as little as a single cell.

Total RNA Purification Plus Micro Kit Contents

1. Buffer RL
2. Wash Solution A
3. Elution Solution A
4. Genomic DNA Removal Column
5. RNA Purification Micro Column
6. Collection Tubes
7. Elution tubes (1.7 mL)
8. Product Insert

Customer-Supplied Reagents and Equipment

- Benchtop microcentrifuge
- 95% ethanol
- β -mercaptoethanol
- PBS (Cultured cells)
- Mortar and pestle, liquid nitrogen, 70% ethanol (Tissue)
- Sterile, fine forceps, 42°C water bath, 70% ethanol (LCM)

Shipping Conditions

The Total RNA Purification Plus Micro Kit is shipped at room temperature.

Storage Conditions

All solutions should be kept tightly sealed and stored at room temperature. These reagents should remain stable for 1 year in their unopened containers.

Cat #	Description	Quantity
48500	Total RNA Purification Plus Micro Kit	50 preps