

VZV TaqMan Probe/Primer and Control Set - 100 Reactions

Product # TM36710

| Component | Product # TM36710 (100 Reactions) |
|------------------------|-----------------------------------|
| VZV Primer & Probe Mix | 1 x 280 µL |
| VZV Positive Control | 1 x 150 µL |
| Nuclease-Free Water | 1 x 1.25 mL |

Applications

TaqMan Probe/Primer and Control Set for VZV detection using real-time PCR based on the use of TaqMan® technology

Storage Conditions and Product Stability

Store all the components at -20°C upon arrival. All kit components can be stored for 2 years after the date of production without showing any reduction in performance.

Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

Customer-Supplied Reagents and Equipment

Appropriate Real-Time PCR Instrument with FAM and HEX (or VIC) filter channel

Procedures

Recommended Reaction Conditions:

- VZV Primer Mix (2 µL per 20 µL PCR reaction)
- VZV Positive Control (2-8 µL per 20 µL PCR reaction)

Description

The positive control is a cloned partial VZV genomic DNA.

Recommended TaqMan PCR Assay

| PCR Cycle | Step | Temperature | Duration |
|----------------------|--------|-------------|----------|
| <i>Cycle 1</i> | Step 1 | 95°C | 3 min |
| | Step 2 | 60°C | 30 sec |
| <i>Cycle 2 (40x)</i> | Step 1 | 95°C | 15 sec |
| | Step 2 | 60°C | 30 sec |

Technical Support

Contact our Technical Support Team between the hours of 8:30 and 5:30 (Eastern Standard Time) at (905) 227-8848 or Toll Free at 1-866-667-4362. Technical support can also be obtained from our website or through email at techsupport@norgenbiotek.com.