




Total RNA Purification

Including microRNA - Without Phenol




Superior microRNA Diversity



*Ask us about **EXTRAClean**

True Total RNA Recovery



- ✓ Isolate total RNA, including microRNA, without using phenol
- ✓ Highly sensitive down to very small inputs
- ✓ RNA, microRNA, DNA and protein isolation from a single sample available
- ✓ Multiple formats available: column, slurry, magnetic bead, 96-well plate.
- ✓ 40+ sample optimized kits to meet all of your RNA and microRNA needs

COLLECTION & PRESERVATION

EXTRACTION

APPLICATION



norgenbiotek.com

Optimized Kits for Any Sample

- ✓ CELLS
- ✓ TISSUE
- ✓ BACTERIA
- ✓ YEAST

- ✓ FUNGI
- ✓ PLANTS
- ✓ WATER
- ✓ SOIL



- ✓ SINGLE CELL
- ✓ EXOSOMES
- ✓ VIRUSES
- ✓ LEUKOCYTE

- ✓ STOOL
- ✓ URINE
- ✓ PLASMA/SERUM
- ✓ FFPE

Animal Tissue RNA Purification Kit

CAT. 25700

Recover all RNA Species from Fibrous and Challenging Tissue Samples



Fatty Tissue RNA Purification Kit

CAT. 36200


- ✓ High quality RNA from high lipid content tissues
- ✓ Brain, adipose, breast, etc.




Cytoplasmic & Nuclear RNA Purification Kit

CAT. 21000

Fractionate Nuclear and Cytoplasmic RNA




RNA Gel



U2 snRNA RT-PCR
(Nuclear-specific)

microRNA Purification Kit

CAT. 21300



Kits Available in Various Formats to Match Any Elution Requirement

Column Format	Elution Range	Use
Mini Spin Column	30-50 μ L	Routine RNA/microRNA Purifications
Micro Spin Column	20-50 μ L	RNA/microRNA Clean-up and Concentration, FFPE RNA
Single Cell Spin Column	8-20 μ L	LCM, Stem Cells/Sorted Cells, CTC, etc.

*Kits also available in 96-well formats for high throughput RNA extractions and automation.

Ordering Information

FFPE RNA and microRNA Extraction and Purification

Kit Name	Size	Cat No.
● Cytoplasmic and Nuclear RNA Purification Kit	50 preps	21000
FFPE RNA/DNA Purification Plus Kit	50 preps	54300
FFPE RNA/DNA Purification Plus Kit Dx	50 preps	Dx54300

Stool RNA Purification

Kit Name	Size	Cat No.
Stool Total RNA Purification Kit	50 preps	49500
Stool Total RNA Purification Kit Dx	50 preps	Dx49500
Stool Nucleic Acid Isolation Kit	50 preps	45600

Ordering Information

Total RNA Purification (including microRNA)

Kit Name	Size	Cat No.
● Total RNA Purification Kit	50 preps	17200
Total RNA Purification Kit Dx	50 preps	Dx17200
● Total RNA Purification Plus Kit	50 preps	48300
Total RNA Purification 96-Well Kit	2 x 96-well	24300
Total RNA Purification Maxi Kit	8 preps	26800
Total RNA Purification Micro Kit	50 preps	35300
Total RNA Purification Plus Micro Kit	50 preps	48500
microRNA Purification Kit	25 preps	21300
Single Cell RNA Purification Kit	50 preps	51800
Animal Tissue RNA Purification Kit	50 preps	25700
Fatty Tissue RNA Purification Kit	25 preps	36200

Magnetic Bead System

Kit Name	Size	Cat No.
Total RNA Purification Kit	50 preps	75400
Total RNA Purification Kit	2 x 96 Well	75500

Blood/Plasma/Serum RNA Purification

Kit Name	Size	Cat No.
● *Plasma/Serum RNA Purification Kit	50 preps	55000
Plasma/Serum RNA/DNA Purification Kit	50 preps	55200
● Plasma/Serum Cell-Free Circulating and Viral Nucleic Acid Purification Kit	50 preps	56300
Leukocyte RNA Purification	50 preps	21200
Leukocyte RNA Purification 96-Well Kit	2 x 96-well	37800
Preserved Blood RNA Purification Kit I (for use with Tempus Blood RNA Tubes)	50 preps	43400
Preserved Blood RNA Purification Kit II (for use with PAXgene Blood RNA Tubes)	50 preps	43500

Isolate RNA, microRNA, DNA and proteins from the Same Sample

Kit Name	Size	Cat No.
● RNA/DNA Purification Kit	50 preps	48700
RNA/DNA Purification Micro Kit	50 preps	50300
RNA/Protein Purification Plus Kit	50 preps	48200
RNA/DNA/Protein Purification Plus Kit	50 preps	47700
RNA/DNA/Protein Purification Plus Micro Kit	50 preps	51600
RNA/DNA/Protein Purification 96-Well Plus Kit	1 x 96-well	51700

Urine RNA and microRNA Purification

Kit Name	Size	Cat No.
● *Urine Cell-Free Circulating RNA Purification Mini Kit	50 preps	56900
● Urine Cell-Free Circulating and Viral Nucleic Acid Purification Mini Kit	50 preps	59900
Urine Exfoliated Cell RNA Purification Kit	50 preps	22550
Urine microRNA Purification Kit	50 preps	29000
Urine Total RNA Purification Maxi Kit (Slurry Format)	variable	29600
Urine Total RNA Purification Maxi Kit Dx (Slurry Format)	variable	Dx29600
Urine Exosome RNA Isolation Kit	50 preps	47200

Plant and Environmental RNA and microRNA Purification

Kit Name	Size	Cat No.
● Plant/Fungi Total RNA Purification Kit	50 preps	25800
Plant/Fungi Total RNA Purification 96-Well Kit	2 x 96-well	31900
Plant microRNA Purification Kit	50 preps	54700
Plant RNA/DNA Purification Kit	50 preps	24400
Soil Total RNA Purification Kit	50 preps	27750
Water RNA/DNA Purification Kit (0.22 µm)	50 preps	26400
Water RNA/DNA Purification Kit (0.45 µm)	50 preps	26450

Cytoplasmic and Nuclear RNA and microRNA Purification


Kit Name	Size	Cat No.
● Cytoplasmic and Nuclear RNA Purification Kit	50 preps	21000

Clean-up and Concentration of microRNA and Total RNA

Kit Name	Size	Cat No.
*RNA Clean-Up and Concentration Kit	50 preps	23600
*RNA Clean-Up and Concentration 96-Well Kit	2 x 96-well	25100
*RNA Clean-Up and Concentration Micro-Elute Kit	50 preps	61000
CleanAll RNA/DNA Clean-Up and Concentration Kit	50 preps	23800
Oligo Clean-Up and Concentration Kit	50 preps	34100

* **EXTRAClean** Version Available

● More Sizes Available



Order Today!

Three ways to order
 Visit norgenbiotech.com
 Call at 1-866-NORGENB or 905-227-8848
 Email at orders@norgenbiotech.com