

## Nuclease-Free Water

Product # 28001, 28002

## Product Insert

### Components

Cat. #	Size
28001	4 x 1.25 mL
28002	100 x 1.25 mL

### Description

Norgen's Nuclease-Free Water is free of all DNases, RNases and nucleic acids and is suitable for all molecular biology applications requiring nuclease-free water including PCR, RT-PCR and real-time PCR. Our water is prepared without the use of DEPC, and is ready to use without further treatment.

### Quality Control

In accordance with Norgen's ISO-certified Quality Management System, this product has passed a quality control assay to verify the absence of detectable levels of RNase, DNase or nucleic acids.

Tech support: [techsupport@norgenbiotek.com](mailto:techsupport@norgenbiotek.com)

Orders: [orders@norgenbiotek.com](mailto:orders@norgenbiotek.com)

Toll Free: 1-866-667-4362 | Phone: 905-227-8848

Fax: 905 227 1061

# Nuclease-Free Water

Catalogue Number: 28001, 28002

## Applications:

For use in any molecular biology application.

## Storage Conditions

Norgen's Nuclease-Free Water should be stored at room temperature. Storage at +4°C or -20°C is optional.

## Precautions and Disclaimers

This product is designed for research purposes only. It is not intended for human or diagnostic use.

ISO 13485:2003, ISO 9001:2008,  
and ISO 15189:2012 Registered

Product of CANADA

NORGEN  
BIOTEK  CORP.